

Baumanometer

ANEROID: Pocket Aneroid 1050

The Baum Pocket Aneroid is a 300 mmHg portable instrument. It's built to agree with our master instruments, which are traceable to the United States National Institute of Standards and Technology. It complies with ANSI/AAMI SP-9 for accuracy and performance (1% full scale +/- 3 mmHg). The unique, high contrast scale and red pointer increase visibility and reduce parallax error. They are manufactured with a precision low friction movement and twin pre-seasoned phosphor-bronze diaphragm capsules.



Guarantee:

Every Baum Aneroid instrument is guaranteed against defects in workmanship and material for a period of one year from date of purchase. This guarantee is invalid when other than genuine W.A. Baum parts are used.

We recommend that every aneroid instrument be checked for accuracy periodically by following our ["Testing Instrument Accuracy"](#) procedure to determine if it needs adjustment. Adjustments will be made to Baum aneroid instruments without charge when returned, postage prepaid unless the instrument has been damaged in shipment, tampered with or abused. Our pledge is that our customers will be completely satisfied with our products.

Ordering Information:

The die cast zinc manometer housing is finished in warm gray textured baked enamel. The standard Adult [Calibrated® V-Lok® Cuff](#) is medium blue Dacron™ polyester. The adult size is also available in dark brown, dark green, scarlet, tan and gray.

Pocket Aneroid

Cat. No.	Description
1050	Pocket Aneroid w/ Adult Calibrated® V-Lok® & case
1050NL	Pocket Aneroid w/ Adult Non-Latex Calibrated® V-Lok® & case
1051	Pocket Aneroid w/ Child/Sm. Adult Calibrated® V-Lok® & case