



## BiPAP<sup>®</sup> Auto M Series with Bi-Flex<sup>®</sup>

BiPAP<sup>®</sup> Auto M Series with Bi-Flex<sup>®</sup>

### Specifications

PRESSURE	4 to 25 cm H <sub>2</sub> O
RAMP TIME	0 to 45 minutes (5-minute increments)
STARTING RAMP PRESSURE	4 cm H <sub>2</sub> O to EPAP; patient-adjustable
DIMENSIONS	7.5"(L) x 5.0"(W) x 3.1"(H) 19.0 cm(L) x 12.7 cm(W) x 7.9 cm(H)
WEIGHT	2.2 lbs. (1.0 kg)
FILTERS	Pollen and optional ultra-fine
DEVICE SET-UP	LCD keypad
DATA STORAGE CAPACITY (minimum)	Display: 7 and 30-day averages SmartCard: 6 months at 3 usage sessions per day, 7 days in-depth analysis
DATA STORAGE CONTENT	SmartCard: Date/time, hours of usage patterns, AHI, leak, snore, FOSQ
COMPLIANCE MONITORING	Breathing detection
ALTITUDE COMPENSATION	Automatic
ELECTRICAL REQUIREMENTS	100 – 240 VAC 50/60 Hz
HUMIDIFICATION	Integrates with M Series Heated and M Series Pass-Over Humidifiers
DC POWER	Direct-connect cord

WARRANTY	2 years
MISC	Lighted LEDs, Auto on/off, mask off alert, patient reminders

## Part Numbers and Accessories

<b>BiPAP Auto SYSTEMS</b>	<b>Part #</b>
BiPAP Auto M Series w/SmartCard	DS700S
BiPAP Auto M Series w/SmartCard (includes humidifier)	DS700HS
<b>M SERIES ACCESSORIES</b>	<b>Part #</b>
POLLEN FILTERS	1029330
ULTRA-FINE FILTERS	1029331
CARRYING CASE	1022219
TRAVEL PACK	1028122
DC ADAPTER CORD	1001956
SHIELDED DC CORD SYSTEM (incl. shielded DC cord and battery adapter cable)	1001979
REUSABLE PATIENT TUBING	622038
AC POWER CORD (5 ft)	1028078
EXTERNAL POWER SUPPLY	1015642
DC Standalone Heated Humidifier	1049109 (compatible with 1003758)
DC Standalone Passover Humidifier	1051081 (compatible with 1003758)
R3B Chamber kit (w/J-Tube openings)	1003758 (compatible with 1049109 and 1037652)
R2B Chamber kit (w/side by side opening)	1003757 (compatible with 1022257)
Universal RP Tan Kit (both tops)	1003756 (compatible with 1049109, 1037652, 1022257)
R3B Inlet/Outlet Seal, Qty. 2	1003784 (compatible with 1003758)
<b>Obsolete Products</b>	<b>Part #</b>
AC Standalone Heated Humidifier	1037652 (compatible with 1003758)
Original M Series Heated Hum.	1022257 (compatible with 1003757)
Original M Series Passover	1029090 (compatible with 1029533)
<b>ENCORE PRO SOFTWARE AND SMARTCARD</b>	<b>Part #</b>
Encore Pro Software	1009583
ENCORE USB SMART CARD READER/WRITER	1014352
Pack of SmartCards (10 cards)	1003304
FOSQ Clinical Reprint (Sleep, 20(10): 835-843)	1000141
FOSQ Questionnaire, English (Call Customer Service for other languages)	1004550