



## BiPAP® Plus M Series with Bi-Flex®

### Specifications

PRESSURE	4 to 25 cm H <sub>2</sub> O
RAMP TIME	0 to 45 minutes (5-minute increments)
STARTING RAMP PRESSURE	4 cm H <sub>2</sub> O to EPAP; patient-adjustable
DIMENSIONS	7.5"(L) x 5.0"(W) x 3.1"(H) 19.0 cm(L) x 12.7 cm(W) x 7.9 cm(H)
WEIGHT	2.2 lbs. (1.0 kg)
FILTERS	Pollen and optional ultra-fine
DEVICE SET-UP	LCD keypad
DATA STORAGE CAPACITY (minimum)	Therapy meter
DATA STORAGE CONTENT	Total therapy hours, Usage sessions > 4 hours
COMPLIANCE MONITORING	Breathing detection
ALTITUDE COMPENSATION	Automatic
ELECTRICAL REQUIREMENTS	100 – 240 VAC 50/60 Hz
HUMIDIFICATION	Integrates with M Series Heated and M Series Pass-Over Humidifiers
DC POWER	Direct-connect cord
WARRANTY	2 years
MISC	Lighted LEDs, Auto on/off, mask off alert, patient reminders

DATA STORAGE CAPACITY (minimum)	Therapy meter
---------------------------------	---------------

DATA STORAGE CONTENT	Total therapy hours, Usage sessions > 4 hours
COMPLIANCE MONITORING	Breathing detection
ALTITUDE COMPENSATION	Automatic
ELECTRICAL REQUIREMENTS	100 – 240 VAC 50/60 Hz
HUMIDIFICATION	Integrates with M Series Heated and M Series Pass-Over Humidifiers
DC POWER	Direct-connect cord
WARRANTY	2 years
MISC	Lighted LEDs, Auto on/off, mask off alert, patient reminders

BiPAP® Plus M Series with Bi-Flex®

## Part Numbers and Accessories

<b>BiPAP Plus SYSTEMS</b>	<b>Part #</b>
BiPAP Plus M Series	DS600
BiPAP Plus M Series (includes humidifier)	DS600H
<b>M SERIES ACCESSORIES</b>	<b>Part #</b>
POLLEN FILTERS	1029330
ULTRA-FINE FILTERS	1029331
CARRYING CASE	1022219
TRAVEL PACK	1028122
DC ADAPTER CORD	1001956
SHIELDED DC CORD SYSTEM (incl. shielded DC cord and battery adapter cable)	1001979
REUSABLE PATIENT TUBING	622038
AC POWER CORD (5 ft )	1028078
EXTERNAL POWER SUPPLY	1015642
DC Standalone Heated Humidifier	1049109 (compatible with 1003758)
DC Standalone Passover Humidifier	1051081 (compatible with 1003758)
R3B Chamber kit (w/J-Tube openings)	1003758 (compatible with 1049109 and 1037652)
R2B Chamber kit (w/side by side opening)	1003757 (compatible with 1022257)
Universal RP Tan Kit (both tops)	1003756 (compatible with 1049109, 1037652, 1022257)
R3B Inlet/Outlet Seal, Qty. 2	1003784 (compatible with 1003758)