



Image shown Aaron 1250 –G OB-G (Available)

Aaron 1250

Available Accessories Include:

- Disposable and Reusable Pencils
see Pencils page
- Bipolar Forceps (Reusable)
see Forceps page
- Disposable Electrodes
see Disposables page
- Reusable Electrodes
see Reusables page
- Disposable Loops
see Loops page
- A1253 Two Pedal Footswitch

- A1255A Adaptor for Connecting Footswitching Pencil (ESP7)
*Sterile
- A1252C Connecting Cord for ESRE and ESRS Grounding Pads
*Sterile
- ESRS Solid Adult Return Electrode (Grounding Pad)
- ESRE Split Adult Return Electrode (Grounding Pad)
- ESRS Solid Adult Return Electrode (Grounding Pad) with 2.8 M Cable, 50/box
- ESRE Split Adult Return Electrode (Grounding Pad) with 2.8 M Cable, 50/box
- ESMS-C Mobile Stand with Instrument Drawer and Bottom Smoke Evacuator Tray
- ESMS Mobile Stand with Instrument Drawer
- ESMS-BT Bottom Smoke Evacuator Tray